

Toshiba Eccocee CX Technical Specifications

Applications

Abdominal
Gynecology
Vascular
Small Parts
Urology

Extended Options

Third Transducer Connector
External Stress Echo
Annular Array Upgrade
Micro and Biopsy Transducers

System Dimensions

33.5" Depth
18.1" Width
51.2" Height
265 lbs. (120kg)

Additional Features

B-Mode and Color Image Progressing
Velocity Range Processing
Auto Focus
Palm Controller
Dual Trackball

Imaging Modes

B-Mode
M-Mode
PW Doppler
Color Flow

Standard Features

12-inch Color Monitor
Color Angio
3D Color Perspective
Digital Image Filter
Triple Frequency Dynamic Beam Control
Real-Time Triplex Mode
Real-Time Pan/Expand
Independent Steering
OB Advisor
Image Memory

Available Transducers

3.0 thru 6.0 MHz Abdominal Convex
3.0 thru 6.0 MHz Abdominal Micro Convex
3.75 MHz Linear (Biopsy)
3.75 thru 5.5 MHz Abdominal Convex
5.0 thru 7.0 MHz Single Plane Intracavity Convex
5.0 thru 7.0 MHz BiPlane Intracavity Convex
5.0 thru 7.0 MHz Endovaginal Convex
5.0 thru 8.0 MHz Vascular Linear
5.0 thru 8.0 MHz Small Parts Linear
6.0 thru 10.0 MHz Small Parts Linear



For more information call 1-866-4-ENCORE